

Right-hand use


Left-hand use

CleverControl Cuboid


Ergonomy and flexibility are key here – control wheels for moving device motors are positioned at ergonomic angles, and the buttons are easy to access.

The button-functionality change depending on the functions displayed on the screen, which gives the user an eight-button user interface, but effectively unlimited functionality from these buttons. The colour display screen shows device coordinates, offers additional device functionality and enables the user to customize device properties. This customization can be linked to user or experiment requirements and saved under different user identities.

Cuboid working screen



PC Software working screen



CleverControl PC Software

The CleverControl PC software offers the same functionality as on the Cuboid software, and more. We have more freedom with space, the positioning and amount of control buttons on the GUI, and also to provide regular updates.